

STUDY SUB-NANOMETER MEMBRANE FLUCTUATIONS IN SINGLE CELLS USING A PLASMONIC IMAGING MICROSCOPE

A THESIS BY
SURAJ KHOCHARE

Committee Members:


- 1) Dr. Xiaonan Shan (Committee Chair), UH ECE
- 2) Dr. Wei-Chuan Shih, UH ECE
- 3) Dr. Ashutosh Agrawal, UH Mech
- 4) Dr. David Mayerich, UH ECE

Outline:

- 1. Background**
- 2. Surface Plasmon Resonance Imaging**
- 3. Sub-Nanometer Membrane Fluctuations**
- 4. Cell Heterogeneity**
- 5. Cell Metastasis**
- 6. Conclusion**
- 7. Future Scope**

Background

World Health Organization Report




According to a survey by American Cancer Society in 2017,

- 15.5 million Americans with a history of cancer were alive on January 1, 2016.
- Around 1.6 million people were expected to be diagnosed with cancer in 2017 (excluding the ones with non-invasive cancer).
- Around 0.5 million were expected to die of cancer in 2017 (1650 people per day).

Do we know the reason for cancer? | If not, can we study cancer cells? | Can we find a remedy?



Introduction to Membrane Fluctuations:

Structure of a Biological Cell



- ✓ Live cells undergo continuous active processes.
- ✓ **Ex. :-** Metabolism, Metastasis, Mitosis, etc.

Cause of Cell Membrane Fluctuations





Active fluctuations:

- ✓ Cross-membrane ion-transport.
- ✓ Cell structure
- ✓ Endocytosis and Exocytosis


Techniques to Image Membrane Fluctuations:

Quantitative Phase Imaging



QPI Image

Interference Reflectance Microscopy



Our Approach: Surface Plasmon Resonance Imaging

Sensitive only to
bottom membrane

- Image only the bottom cell membrane.
- High sensitivity in z-direction, less distortion.


High Sensitivity :
Sub-nanometer
Membrane Fluctuations

- Study metabolism.
- Study metastasis.

Single Cells

- Study fundamental properties.
- Study cell heterogeneity.


Principle of Surface Plasmon Resonance (SPR) Imaging



SPR Interface

$$\sin \theta_{sp} = \frac{\sqrt{\frac{\epsilon_m n_s^2}{\epsilon_m + n_s^2}}}{n_g}$$



SPR Angle




SPR Response

✓ High Sensitivity

Live Cell Imaging Using SPR






Transmitted Image



SPR Image



Mapping Reflected Light Intensity to Distance from Substrate



Calibration Curve

Simulation Model

Layer 5

Cytoplasm

Layer 4

Cell Membrane

Layer 3


Culture Medium

Layer 2




Au

Layer 1


BK7



Imaging the Sub-nanometer Membrane Fluctuations



Analyse Sub-Nanometer Membrane Fluctuations





- ✓ Image small yet swift movement.
- ✓ Distribution range for fluctuations.
- ✓ These movements may contribute to respiration, metabolism, active fluctuations and thermal fluctuations.

Confirm Active Membrane Movement By Cell Fixation

- Add Paraformaldehyde (4% in PBS)









Add 100 μ l PFA in 200 μ l Culture Medium (50% Conc.)



Record Entire Fixation Process


- Membrane Fluctuations: Brownian or Active?
- Highly concentrated dose kills the cells.
- Therefore, expect no membrane fluctuations.
- This confirms that we image active movement.

Quick Review




Cell Heterogeneity


a)



b)






c)




- ✓ Characteristics of a population may not reflect characteristics of an individual cell.
- ✓ Different cells are at different stages of growth at a given time point.
- ✓ **Factors Involved:** Environmental Conditions, Genetic Variations, interaction with the surrounding.

Intra-Cell Correlation and Differences

HeLa Cells




HT-1080 Cells







- ✓ Center is more active than edges.
- ✓ Observe energy expenditure.

Intra-Cell Transient Heterogeneity

HeLa Cells





- Analyse sub-regions ($5 \times 5 \mu\text{m}$).
- Transient heterogeneity.
- Related to active processes.



Long Duration Cell Membrane Movement


(Scale Bar = 10 μm)

3D Video Displaying Membrane Movement




- ✓ We see a collective movement (10's of nanometers) as we record for longer time.
- ✓ This can be related to physiological processes like metastasis.
- ✓ Do they have a correlation with the sub-nanometer membrane fluctuations?

Relation Between Membrane Movement and Sub-nanometer Fluctuations




- Both movements are correlated to some extent.


Cell Metastasis



HT-1080




Wound Healing Assay




- ✓ HT-1080: Human Fibrosarcoma Cell Line.
- ✓ Highly active cell membrane fluctuations.

Membrane Movement in Metastatic Cells



Bright-Field Image



SPR Image



Membrane Movement




Displaced Area \propto Membrane Movement

Conclusion:

- ✓ Live cell imaging of cell bottom membrane with high resolution in z-direction and less noise.
- ✓ Imaging the cell membrane fluctuations important to study metabolism.
- ✓ We observe heterogeneity in cells which corresponds to multiple factors.
- ✓ Enable to analyse metabolic response of a cell during anti-cancer drug treatments.
- ✓ Cell metastasis can be evaluated based on membrane fluctuations to some extent.

Future Plan

- 1) Using membrane fluctuations as a characteristics to study metabolic activities at single cell level.



- 2) Analyse response of cells to anti-cancer drugs targeting active processes.

Thank You...!!!